

# eyeSeat/Seat assistant

The intelligent seat-analysis system




eyeSeat/Seat assistant is an optional feature available on the @Just swivel-chair model (AJ 26820). Four pressure sensors in the seat and two in the backrest show **how the user is sitting on the chair and the postures they adopt**. For example, if the test person spends a relatively long period of time sitting on the front edge of the seat or without making contact with the backrest, a message reminding them of the need to correct their posture is shown on their PC monitor. An electronic operating manual rounds off this option.



The eyeSeat /Seat assistant registers and documents the test person's seated postures. Behavioural patterns are shown on the monitor and suggestions for improvements are given where necessary, ultimately with a view to achieving a change in behaviour.

The key objective of this system is to ensure that the user achieves a balanced seated posture by adjusting the chair correctly, making full use of its features and alternating between movement and rest.



Ergonomics

## Ergo-Contact

➔ **Malte Lenkeit**/Ergonomics Consultant

Telephone: +49 (0) 30/43 55 76-0, Mobile: +49 (0) 172/869 12 53, E-mail: malte.lenkeit@dauphin.de

➔ **Matthias Thumshirn**/International Head of Training

Telephone: +49 (0) 91 58/17-391, Mobile: +49 (0) 1 72/869 12 84, Fax: +49 (0) 91 58/17-780, E-mail: matthias.thumshirn@dauphin.de



**DAUPHIN**

HumanDesign® Company